

Mycoplasma Removal Reagent (50x) LTAAnB-L03

Product Description

Mycoplasma, first discovered in 1898, was firstly classified as virus. There are more than 100 known species belonging to the genus Mycoplasma. They are the smallest free-living, self-replicating organisms and lack a cell wall so that antibiotics such as penicillin, which interfere with murein formation, are not effective. If not treated, mycoplasma contaminations can lead to diminished cell growth and to the loss of cultures.

MycroXpert is a high efficient antibiotic active at low concentration for a broad range of mycoplasma subspecies. It acts by interference with the DNA replication machinery and thus inhibiting mycoplasma growth. **Mycoplasma Removal Reagent** is provided as a 50x concentrate.

Product Specifications

pH	4.0 – 4.5
Endotoxin	≤10 EU/ml
Cell Culture	Tested
Sterility	Tested
Storage	Store at ≤-15°C. Avoid repeated freeze thaw cycles.

Procedure

Treatment with **Mycoplasma Removal Reagent** includes three or more cycles (passages). It is recommended to check the culture every week for presence of mycoplasma (e.g. ELISA, fluorescence staining or PCR). Please note, that using PCR for detection may lead to false positive results during treatment with **Mycoplasma Removal Reagent** since this method cannot distinguish between living and dead mycoplasma. We recommend using PCR as detection method only approx. 2 weeks after treatment with **Mycoplasma Removal Reagent** in order to avoid false positive results.

Normally mycoplasma will be removed after 2 to 3 weeks.


Working concentration

Mycoplasma Removal Reagent is provided as a 50x concentrate.

- Working concentration: 1x.
Therefore dilute the stock solution 1 : 50 (Dilution 1)
Example: Add 10 ml **Mycoplasma Removal Reagent** to 500 ml base medium.
- If required increase **Mycoplasma Removal Reagent** concentrations as follows in different subsequent steps:
Dilution 2: 1 : 40
Dilution 3: 1 : 30
Dilution 4: 1 : 25

Instruction for Use

- Dilution 1 (1 : 50)
- Remove culture medium from culture vessels, wash cells with PBS and detach cells using Trypsin or Accutase.
- Count cells and plate in fresh medium (with **Mycoplasma Removal Reagent**).
- Cultivate cells for 2 to 3 days in the usual way.



UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918,
V.Į. Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27

- Repeat the passaging cycle at least twice.
- Check for mycoplasma contamination (Staining method, cultivation or ELISA).
- If mycoplasma positive, check methods (positive and negative controls and evaluate).
- If required repeat treatment with increased **Mycoplasma Removal Reagent** concentrations (Dilution 2, Dilution 3 or Dilution 4).

Precautions and Disclaimer

- This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.