



UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registru centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Bovine Serum Albumin, Fatty Acid Free, Origin USA 100g LTSAAF100 Bovine

Serum Albumin, Fatty Acid Free, Origin USA 50g LTSAAF50

General Information

Bovine Serum Albumin (also known as BSA) is a serum albumin protein derived from bovine blood. It is often used as a protein concentration standard in lab experiments.

BSA Fatty Acid Free has extremely low levels of total lipids and fatty acids and makes it suitable for various applications in which lipids and fatty acids may interfere.

BSA has numerous biochemical applications including ELISAs (Enzyme-Linked Immunosorbent Assay), immunoblots, and immunohistochemistry. It is also used as a stabilizer and cell growth supplement in cell and microbial cultures.

Applications

- Additive for cell culture
- Binding and transport studies
- Protein stabilizer
- Blocking agent

Appearance	White or off-white powder
Application	Before Application in cell culture, prepare a sterile filtered stock solution in PBS or water.
CAS No.	9048-46-8
Molecular weight	66 kDa
Storage	Store at +2°C to +25°C in dry conditions protected from light.
Shipping conditions	Ambient

Precautions and Disclaimer

This product is for research use only.