



## UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,  
Rugijų 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.  
Registrų centras Vilniaus filialas  
Tel/fax +370 5 216 02 27

### Dulbecco's PBS(10x), w/o Ca & Mg, w/o Phenol Red LT BSS05

#### General Information

Dulbecco's Phosphate Buffered Saline (DPBS) is a water-based salt solution used in various research areas. The buffer is non-toxic and provides cells with water and inorganic ions which are necessary for cell metabolism. It can be also used for rinsing, as a diluent or reference solution.

#### Product Specifications

Appearance	Clear liquid
Storage and shelf life	Store at +2 to +25°C upon arrival. Once opened store at +2 to +25°C and use within 24 months.
Shipping conditions	Ambient

#### Formulation

Component	concentration mg/l
Potassium chloride	200
Potassium phosphate monobasic anhydrous	200
Sodium chloride	8000
Sodium phosphate dibasic anhydrous	1150

#### Precautions and Disclaimer

This product is for research use only.