



## UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,  
Rugijų 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.  
Registrų centras Vilniaus filialas  
Tel/fax +370 5 216 02 27

### Hanks'BalancedSalts,w/oCa&Mg,withPhenolRed,sterile- filtered LTBS03

#### General Information

Hank'sBalancedSaltSolution(HBSS)isusedforavarietyofcellcultureapplications,suchas rinsingcells,transportingcellsortissuesamples, dilutingcellsforcounting,andpreparing reagents.TheessentialfunctionofabalancedsaltsolutionistomaintainpHandosmoticbalanceas well as providing cells with water and essential inorganic ions.

#### Product Specifications

Appearance	Clear red liquid
Storage and shelf life	Storeat+2to+8°Cuponarrival. Onceopenedstoreat+2to+8°Candusewithin24month.
Shipping conditions	Ambient

#### Formulation

Component	concentration mg/l
D-Glucose	1000
KCl	400
KH <sub>2</sub> PO <sub>4</sub>	60
NaCl	8000
NaHCO <sub>3</sub>	350
Na <sub>2</sub> HPO <sub>4</sub>	48
Phenol red sodium salt	11

#### Precautions and Disclaimer

Thisproductisforresearchuseonly.