



UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,
Rugijų 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27

DMEMHighGlucose(4.5g/l),withL-GlutamineLTLM-B03

General Information

Appearance	Clear red orange solution
CO ₂ concentration,optimum	8.5%
Storageandshelflife	Storeat+2°to+8°Cprotectedfromlight. Onceopened,storeat4°Candusewithin6-8weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/l	Concentration mg/l
AminoAcids		D-Calcium Pantothenate 4.00
Glycine	30.00	FolicAcid 4.00
L-Arginine HCl	84.00	myo-Inositol 7.20
L-Cystine 2 HCl	48.00	Nicotinamide 4.00
L-Glutamine	584.00	Pyridoxal HCl 4.00
L-HistidineHClH ₂ O	42.00	Riboflavin 0.40
L-Isoleucine	105.00	ThiamineHCl 4.00
L-Leucine	105.00	Inorganic Salts
L-Lysine HCl	146.00	CaCl ₂ 2H ₂ O 265.00
L-Methionine	30.00	Fe(NO ₃) ₃ 9H ₂ O 0.10
L-Phenylalanine	66.00	KCl 400.00
L-Serine	42.00	MgSO ₄ 97.67
L-Threonine	95.00	NaCl 6400.00
L-Tryptophan	16.00	NaHCO ₃ 3700.00
L-Tyrosine 2 Na	89.47	NaH ₂ PO ₄ 109.00
L-Valine	94.00	Other Components
Vitamins		D-Glucose 4500.00
Choline chloride	4.00	Phenol Red Sodium Salt 15.90

Precautions and Disclaimer

Thisproductisforresearchuseonly.