



UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,
Rugijų 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27

DMEM High Glucose (4.5g/l), with Stable Glutamine, with 25mM HEPES LTLM-B04

General Information

Appearance	Clear red orange solution
CO ₂ concentration, optimum	8.5%
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened, store at 4°C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/l	Concentration mg/l
Amino Acids		D-Calcium Pantothenate 4.00
Glycine	30.00	Folic Acid 4.00
L-Alanyl-L-Glutamine	862.00	myo-Inositol 7.20
L-Arginine HCl	84.00	Nicotinamide 4.00
L-Cystine 2 HCl	48.00	Pyridoxal HCl 4.00
L-Histidine HCl·H ₂ O	42.00	Riboflavin 0.40
L-Isoleucine	105.00	Thiamine HCl 4.00
L-Leucine	105.00	Inorganic Salts
L-Lysine HCl	146.00	CaCl ₂ ·2H ₂ O 265.00
L-Methionine	30.00	Fe(NO ₃) ₃ ·9H ₂ O 0.10
L-Phenylalanine	66.00	KCl 400.00
L-Serine	42.00	MgSO ₄ 97.67
L-Threonine	95.00	NaCl 6400.00
L-Tryptophan	16.00	NaHCO ₃ 3700.00
L-Tyrosine 2 Na	89.47	NaH ₂ PO ₄ 109.00
L-Valine	94.00	Other Components
Vitamins		D-Glucose 4500.00
Choline chloride	4.00	HEPES 5958.00
		Phenol Red Sodium Salt 15.90

Precautions and Disclaimer

This product is for research use only.