



## UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,  
Rugijų 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.  
Registrų centras Vilniaus filialas  
Tel/fax +370 5 216 02 27

### DMEMLowGlucose(1g/l),w/oL-Glutamine,withSodiumPyruvate LTLM-B12

#### General Information

Appearance	Clear red orange solution
CO <sub>2</sub> concentration,optimum	8.5%
Storageandshelflife	Storeat+2°to+8°Cprotectedfromlight. Onceopened,storeat4°C andusewithin6-8weeks.
Shipping conditions	Ambient

#### Formulation

Components	Concentration mg/l	Concentration mg/l
<b>AminoAcids</b>		D-Calcium Pantothenate 4.00
Glycine	30.00	FolicAcid 4.00
L-Arginine HCl	84.00	myo-Inositol 7.20
L-Cystine 2 HCl	48.00	Nicotinamide 4.00
L-HistidineHClH <sub>2</sub> O	42.00	Pyridoxal HCl 4.00
L-Isoleucine	105.00	Riboflavin 0.40
L-Leucine	105.00	ThiamineHCl 4.00
L-Lysine HCl	146.00	<b>Inorganic Salts</b>
L-Methionine	30.00	CaCl <sub>2</sub> 2H <sub>2</sub> O 265.00
L-Phenylalanine	66.00	Fe(NO <sub>3</sub> ) <sub>3</sub> 9H <sub>2</sub> O 0.10
L-Serine	42.00	KCl 400.00
L-Threonine	95.00	MgSO <sub>4</sub> 97.67
L-Tryptophan	16.00	NaCl 6400.00
L-Tyrosine 2 Na	89.47	NaHCO <sub>3</sub> 3700.00
L-Valine	94.00	NaH <sub>2</sub> PO <sub>4</sub> 109.00
<b>Vitamins</b>		<b>Other Components</b>
Choline chloride	4.00	D-Glucose 1000.00
		Phenol Red Sodium Salt 15.90
		Sodium Pyruvate 110.00

#### Precautions and Disclaimer

Thisproductisforresearchuseonly.