

UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registru centras Vilniaus filialas
Tel/fax +370 5 216 02 27

MEMwithEarle'sSalts,withL-GlutamineLTLM-B18


General Information

Appearance	Clearredorangesolution
CO ₂ concentration, optimum	5.0%
Storage and shelf life	Storeat+2°C–+8°Cprotectedfromlight. Onceopened,storeat4°Candusewithin6-8weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/l
Amino Acids	
L-Arginine HCl	126.00
L-Cystine 2 HCl	31.30
L-Glutamine	292.00
L-Histidine HCl H ₂ O	42.00
L-Isoleucine	52.00
L-Leucine	52.00
L-Lysine HCl	72.50
L-Methionine	15.00
L-Phenylalanine	32.00
L-Threonine	48.00
L-Tryptophan	10.00
L-Tyrosine 2 Na 2 H ₂ O	51.90
L-Valine	46.00
Vitamins	
Choline chloride	1.00
D-Calcium Pantothenate	1.00
Folic Acid	1.00
myo-Inositol	2.00
Nicotinamide	1.00
Pyridoxal HCl	1.00
Riboflavin	0.10
Thiamine HCl	1.00
Inorganic Salts	
CaCl ₂ 2 H ₂ O	265.00
KCl	400.00
MgSO ₄	97.67
NaCl	6800.00
NaHCO ₃	2200.00

NaH ₂ PO ₄	122.00
Other Components	
D-Glucose	1000.00
Phenol Red Sodium Salt	11.00



UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registru centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Precautions and Disclaimer

This product is for research use only.